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Compass™ CR Series Scale Quick Start Guide

Please visit www.ohaus.com/support to download user manual.

SAFETY PRECAUTIONS

- Only operate the scale in dry locations.
- Do not operate the scale in hazardous locations.
- Do not drop or overload the scale.
- When operating the scale from mains power, only use the AC Adapter specified by Ohaus.

INSTALLATION

- Remove the plastic transportation protective insert (0.1g model only) for use. Save the plastic insert to ensure safe transport. See figure 1.
- Remove the battery cover on bottom of the scale and insert the 3 LR6 (AA) batteries into the battery compartment. See figure 2.

TWO BUTTON KEYPAD

- Press  to turn the scale on.
- Long press  to turn the scale off.
- Press  to set zero.
- Press  to change units.
- Long press  to start calibration.

CALIBRATION

1. Long press  until [CR L] is displayed.
2. Press  to set the zero point.
3. When [Cxxx] is displayed, place xxx g weight on the platform.
4. Press  to set the span point.
5. After calibration, the display returns to weighing.

Note:

[CR L E] will appear if the calibration steps are not followed or the wrong weight was used.

Turning the scale off aborts calibration.

MAINTENANCE

- The exterior surfaces of the instrument may be cleaned with a cloth dampened with water and a mild detergent. Do not use solvents, harsh chemicals, ammonia or abrasive cleaning agents.
- Register your scale under *Warranty Registration* at www.ohaus.com/support.

For technical issues contact an Authorized OHAUS Service Agent. Please visit our website www.ohaus.com to locate the OHAUS office nearest you.

COMPLIANCE

	<p>This product complies with Directives 2011/65/EU (RoHS), 2014/30/EU (EMC) and 2014/35/EU (LVD). The EU Declaration of Conformity is available online at www.ohaus.com/ce.</p>
	<p>This product complies with the EU Directive 2012/19/EU (WEEE). For disposal instructions in Europe, refer to www.ohaus.com/weee</p>

FCC NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help.
- Changes or modifications not expressly approved by OHAUS Corporation could void the user's authority to operate the equipment.

Figure 1 / Figura 1 / Figure 1



Figure 2 / Figura 2 / Figure 2

**Technical Data / Datos Técnicos / Données Techniques**

Models	Max	d
CR221	220 g	0.1 g
CR621	620 g	0.1 g
CR2200	2200 g	1 g
CR5200	5200 g	1 g

Ohaus Corporation
7 Campus Drive
Suite 310
Parsippany, NJ 07054 USA
Tel: +1 973 377 9000 ·
Fax: +1 973 944 7177

www.ohaus.com

* 3 0 4 6 7 7 4 9 *

Please scan the below QR code to download the user manual: